

The valves perform within a tolerance range of the average pressure as shown here regardless of gravitational orientation:

Settings 1, 2, 3 ± 20 mmH₂O



Setting 4 ± 25 mmH₂O

Settings 5, 6, 7 ± 35 mmH₂O

Setting 8 is intended to limit flow through the valve and has a minimum pressure of 400 mmH₂O.

Codman® CERTAS™ Plus Programmable Valves were made available after March 2015

X-Ray Performance Setting Guide

-  Right Hand Side (RHS) X-Ray Marker
-  Setting Indicator

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

• Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

• Warning: Applicable laws restrict these products to sale by or on the order of a physician.

• Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 • 888-980-7742 fax

International +1 609-936-5400 • +1 609-750-4259 fax

integralife.com

INTEGRA
LIMIT UNCERTAINTY

Codman, Integra, and the Integra Logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Certas is a trademark of Integra LifeSciences Corporation or its subsidiaries. ©2018 Integra LifeSciences Corporation. All rights reserved. Printed in USA. OM 0880331-2-EN